

FastStart® Electrodes

Interferential and Neuromuscular Stimulation Electrodes

FastStart electrodes are among the most efficient electrodes on the market.

- FastStart electrodes are specifically designed to maximize single patient reuse
- Silver/Silver Chloride ensures optimal current distribution over the entire surface of the pad
- Lift tab provides easy, non-contact application and removal
- New backing substrate conforms evenly to body contours
- Can be effectively used with heat or cold therapy

FastStartSM
Orthopedic Products

800.652.1135

FastStart® Electrodes

Product Profile

Ordering Information



Color-coded for
interferential stimulation



Unique backing conforms
to body contours

Product Number	Description
750-2020	2.0" X 2.0" Electrodes 20 Pouches of 4 each
750-1535	1.5" X 3.5" Electrodes 20 Pouches of 4 each
750-3050	3.0" X 5.0" Electrode 20 Pouches of 2 each
750-0030	3.0" Round Electrodes 20 Pouches of 4 each
750-0020	2.0" Round Electrodes 20 Pouches of 4 each
750-1515	2.2" X 2.5" Electrodes Pouch of 40 each

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of the device.

A Self-adhesive Single Patient Reusable Neuromuscular Stimulation Electrode

- FastStart electrodes are specifically designed to maximize single patient reuse in the harshest application environment
- FastStart electrodes use a Silver/Silver Chloride coated conductive substrate to provide uniform distribution of electrical current over the entire surface of the electrode for optimal stimulation targeted muscle groups
- A laminated cloth backing conforms to body contours, especially around areas with a drastic radius. The FastStart electrode remains in place during all forms of rehabilitation therapy.
- FastStart electrodes are extremely stable, regardless of the degree of patient activity
- The electrode can also be safely and effectively used during hot and cold therapy
- With six sizes to choose from, the right size exists for every patient
- A lift tab assists in removal of the electrode
- The resealable pouch provides easy storage for same patient reuse
- Perspiration and water enhance the adhesion of the electrode during high intensity therapy

18022 Cowan Street, Suite 200F, Irvine, CA 92614
800.652.1135 Fax 800.652.1126 www.faststartortho.com

FastStart is a registered trademark and service mark of VisionQuest Industries, Inc. © 2004

FastStartSM
Orthopedic Products

FS04025200