

FreeStride OA™

Custom Osteoarthritis Knee Brace

You've decided to prescribe knee bracing for unicompartmental osteoarthritis. The most important decision has to be made by your patients – they must choose to wear the brace.



Relief Begins with Compliance

While some osteoarthritis patients receive immediate pain relief from bracing, many must wear a brace several weeks to allow inflammation to subside. Early compliance is a crucial step in OA bracing therapy.

The custom FreeStride OA™ brace represents the state of the art in low-profile, comfortable bracing. The new semi-circumferential shells reduce bulk by more than 50%. In fact, the FreeStride OA can be folded to store in a thin briefcase or purse when not in use. This revolutionary design also eliminates the tibial crest pain created by other braces. Soft suede liners extend beyond the shell borders to comfortably keep the brace in place over many hours of continual use.

Patented Pain Reduction

The pain-reducing action of the FreeStride OA brace is maximized at the terminal extension phase of gait, when knee surfaces are under load and most susceptible to wear. This patented brace design ensures that the tibial lever arm is positioned down the mid-medial or mid-lateral line of the tibia during the mid-phase of gait, providing optimal corrective forces.



FreeStride OA™

Available in 12 custom colors



Patient Fit Guarantee

The FreeStride OA™ brace is guaranteed to fit the patient for 90 days following patient acceptance.

To obtain prescription forms or order a FreeStride OA brace, please contact your dedicated VQ OrthoCare representative or call 800.652.1135.

Ease of Use

First-time brace wearers have little patience for complicated fitting procedures. Obese patients may lack the flexibility to thread D-rings. Geriatric patients may not have the strength or dexterity to put on a large, heavy brace.

To address these needs, we started with an ultra-lightweight design that's easy to handle. Elasticized upper and lower straps allow the patient to quickly secure the brace in place, while the rigid mid-straps maximize corrective forces at the knee. Side buckles with push-button release clips make it a snap to put on or remove the brace.

Features

- Custom-contoured hinge profile for accurate application of corrective forces
- Custom formed thigh and tibia shells for optimum fit and corrective forces
- Custom applied corrective forces
- Ex-Clips provide 5, 10, 15, 20 or 25 degree extension stops
- Custom applied tibial rotation to maximize force at terminal extension
- Custom strap lengths for accurate patient fit and ease of application
- Ultra lightweight (13 oz. short medium; 15 oz. standard medium)
- High-strength 6061-T5 aluminum hinge with built-in tibial rotation mechanism
- Fiber-reinforced composite shells are strong and lightweight
- Elasticized hook and loop straps (upper thigh and lower tibia) aid suspension throughout gait
- Non-elasticized mid-straps maximize corrective forces
- Side buckles ease brace application and removal
- Soft suede liners provide comfort and suspension
- Optional suspension wraps to contain excess soft tissue
- Brace can be made using VQ OrthoCare™ Digi-Cast™ digital measurement system or plaster

Ordering Information

	Left	Right
Brace	50607-00-10	50607-00-20
Accessories		
Suspension Wrap, Small (20")	50714-00-02	
Suspension Wrap, Medium (24")	50714-00-03	
Suspension Wrap, Large (30")	50714-00-04	
Suspension Wrap, X-Large (36")	50714-00-05	
Undersleeve, XS	50731-00-01	
Undersleeve, Small	50731-00-02	
Undersleeve, Medium	50731-00-03	
Undersleeve, Large	50731-00-04	
Undersleeve, X-Large	50731-00-05	
Undersleeve, XXL	50731-00-06	
Undersleeve, XXXL	50731-00-07	
Oversleeve	07-50728-00	

Indications

- Mild to severe unicompartmental osteoarthritis

Warnings

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult a physician or knee brace provider.

Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

CAUTION

Federal law restricts this device to sale by, or on the order of, a practitioner licensed by the law of the State in which he/she practices to use or order the use of this device.