

## Prefabricated Ligament Knee Brace

Whether you brace post-surgically, conservatively or for contact sports, your patients need a knee brace that protects their ligaments. They want a brace that is comfortable, doesn't slide down their leg, and doesn't impede their motion – even during full activity.



The patented Active Thigh Cuff (ATC) is less restrictive than rigid thigh braces, allowing full flexion of the leg.

## A Breakthrough in Suspension

The new prefabricated Hawk™ ligament brace combines the benefits of rigid bracing protection with a new level of suspension and comfort.

Featuring the patented Active Thigh Cuff (ATC) system, the Hawk absorbs soft-tissue profile changes during activity. This unique feature allows the brace to exert a consistent skeletal force at the thigh, providing unparalleled suspension and comfort. The strong and lightweight fiber-reinforced composite shells and high-strength 6061-T5 aluminum hinge provide the strength and protection needed by patients with damaged, deficient or repaired ligaments.

## Active Thigh Cuff (ATC) System



These photos demonstrate the action of the ATC system under various conditions (the active straps are in red)



To obtain prescription forms or order a Hawk brace, please contact your dedicated VQ OrthoCare representative or call 800.652.1135.



The new Hawk® ligament brace features the patented Active Thigh Cuff (ATC) system to dynamically suspend the brace from full flexion to full extension.

### Features

- Patented Active Thigh Cuff (ATC)
- Strong and lightweight fiber-reinforced composite shells
- High-strength 6061-T5 aluminum hinge
- Choice of extension stops in 5, 10, 15, 20 or 25 degrees, or full ROM
- Soft suede liners for comfort and suspension
- Selection of right or left, and standard or short length
- Available in Ultra Black finish

### Ordering Information

Hawk – Standard Length	Left	Right
Small	52405-00-12	52405-00-22
Medium	52405-00-13	52405-00-23
Large	52405-00-14	52405-00-24
X-Large	52405-00-15	52405-00-25
XXL	52405-00-16	52405-00-26
Hawk – Short Length	Left	Right
Small	52410-00-12	52410-00-22
Medium	52410-00-13	52410-00-23
Large	52410-00-14	52410-00-24
X-Large	52410-00-15	52410-00-25
XXL	52410-00-16	52410-00-26
Accessories		
Combined Instabilities Strap	58071-00-00	
Undersleeve, XS	50731-00-01	
Undersleeve, Small	50731-00-02	
Undersleeve, Medium	50731-00-03	
Undersleeve, Large	50731-00-04	
Undersleeve, X-Large	50731-00-05	
Undersleeve, XXL	50731-00-06	
Undersleeve, XXXL	50731-00-07	
Oversleeve	07-50728-00	
Patella Guard, Small/Medium	58003-00-03	
Patella Guard, Large/X-Large	58003-00-04	

### Sizing

Size	Thigh Circumference <sup>1</sup>	Calf Circumference <sup>2</sup>
Small	15 1/2" - 18"	12" - 13 3/4"
Medium	18" - 22"	13 3/4" - 15 1/2"
Large	22" - 24 1/2"	15" - 16 1/2"
X Large	24 1/2" - 27 1/2"	16 1/2" - 19 1/4"
XX Large	27 1/2" - 30 1/2"	19 1/4" - 24 1/4"

<sup>1</sup> Measured 6" above mid-patella; <sup>2</sup> Measured 6" below mid-patella

### Warnings

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult a physician or knee brace provider.

### Indications

- ACL, MCL, PCL, LCL and combined instabilities
- Pre- and post-surgery
- Conservative treatment of ligament damage or deficiency
- Prophylactic bracing for contact sports

### Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort, there is no guarantee that injury will be prevented through the use of this product.

#### CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by, or on the order of, a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.