

Ankle Supports

Sport 8 Ankle Stabilizer



07-20001-XX

The Sport 8 Ankle Stabilizer is a convenient and effective stabilizer for all your training needs...offering unprecedented ease of application and support. "Figure 8" straps start and stop at the ankle joint, providing constant resistance to inversion without hindering ankle movement during activity.

Universal Ankle Brace



3-81330

Ideal for acute and chronic ankle instabilities, as well as prophylactic use. Rigid thermoplastic one-piece outer shell helps prevent inversion/eversion. Cushioned liner incorporates neoprene forefoot closure. Elastic proximal strap provides additional compression as needed. Universal design fits right or left.

Ankle Supports

Sport 8 Ankle Stabilizer

Features and Benefits

- Figure 8 straps provide constant resistance to inversion without hindering movement during activity
- Lace-up design for optimal fit and adjustability
- Supra-malleolar strap locks the isometric stirrups for added support
- Patent #5899872

Product #	Size	Measurement*
07-20001-00	XX-Small	10"
07-20001-01	X-Small	10" - 11"
07-20001-02	Small	11" - 12"
07-20001-03	Medium	12" - 13"
07-20001-04	Large	13" - 14"
07-20001-05	X-Large	14" - 15"
07-20001-06	XX-Large	15" - 16"

*Measure circumference of ankle joint at widest point (around heel)

Universal Ankle Brace

Features and Benefits

- Rigid thermoplastic one-piece outer shell helps prevent inversion/eversion
- Cushioned liner incorporates neoprene forefoot closure
- Elastic proximal strap provides additional compression
- Universal design fits right or left ankle
- Patent #5501659

Product #	Size
3-81330	U

Precautions

These products are to be fitted initially by a healthcare professional who is familiar with the purpose for which they are used. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare professionals, care providers involved in the patient's care, and the patient. If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, consult a physician.

Indications

Chronic Ankle Instability, Acute Ankle Sprains, Prophylactic Ankle Support.

Warnings

Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene byproducts.

Notice

While every effort has been made in state-of-the-art techniques to obtain the maximum compatibility of function, strength, durability and comfort. There is no guarantee that injury will be prevented through the use of this product.

CAUTION

Applicable federal, state, and/or local regulations may restrict this device to sale by, or on the order of, a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.