



Functional Arm Brace



A complete upper extremity bracing solution.

- Quick adjust hand piece allows for supinated, neutral and pronated hand positions
- Lightweight bracing alternative to casting
- Range of motion and fixed position at 30 degree increments
- Allows early passive mobilization and rapid progression to active mobilization
- Use in the clinical management of:
 - Humeral & forearm fractures
 - Postoperative shoulder management
 - Biceps reattachment
 - Ulnar nerve transposition

FastStartSM
Orthopedic Products

800.652.1135



Functional Arm Brace

Product Profile



Ordering Information

- The humeral brace can be easily attached to the standard brace on site. Purchase both a standard brace and a humeral brace for on-site modification.
- The humeral brace is designed for use with compression stockinette. Liners can be supplied to special order only.
- The immobilization strap allows complete shoulder immobilization.

Sizes

Forearm circumference is used as a guide for the standard upper limb brace. Humeral circumference is used as a guide for the humeral brace and humeral brace attachment.

	Forearm Circ.	Distal Humeral Circ.	Prox Humeral Circ.
Small	10-13" (25-33 cm)	<9" (23 cm)	<11" (28 cm)
Medium	12-15" (30-38 cm)	9-11" (23-28 cm)	11-13" (28-33 cm)
Large	14-17" (35-43 cm)	>11" (28 cm)	>13" (33 cm)
X-Large	16-19" (40-48 cm)		

Codes

	Small	Medium	Large
Upper limb brace - fixed/neutral	07-30300-00-02	07-30300-00-03	07-30300-00-04
Supination/pronation version	07-30312-00-02	07-30312-00-03	07-30312-00-04
Humeral brace/attachment	07-30301-00-02	07-30301-00-03	07-30301-00-04
Humeral brace pads	07-30309-00-02	07-30309-00-03	07-30309-00-04
Anterior fracture shells (universal size)	07-30302-00-00		

18022 Cowan Street, Suite 200F, Irvine, CA 92614
 800.652.1135 Fax 800.652.1126 www.faststartortho.com

EXO is a trademark of VisionQuest Industries, Inc.
 FastStart is a registered trademark and service mark of VisionQuest Industries, Inc. © 2004

FastStartSM
 Orthopedic Products