



# Functional Knee Brace



For the management of ACL, PCL,  
MCL, LCL and combined instabilities.

- Precise polycentric hinge
- Contoured straps
- Tibial lock
- Strong, lightweight composite construction

**FastStart**<sup>SM</sup>  
Orthopedic Products

800.652.1135




# Functional Knee Brace

## Product Profile

### Features and Benefits

- Our precise polycentric hinge is easy to adjust to the desired range of motion
- Extension stops at 0, 5, 15, 30, 60 and 90 degrees
- Flexion stops at 0, 15, 30, 60 and 90 degrees
- Convenient snap-on hinge covers for easy access
- Contoured straps for longer wear and convenience
- Unique tibial lock for brace stability and constant tibial force
- Standard color is black



Designed to accommodate a variety of sizing requirements, the brace can be utilized by men and women from ages 12+.

#### Size Guidelines

The key measurement is the thigh circumference six inches above the mid patella.

	Inches	Centimeters
Small	15.5-18"	38-44
Medium	18-22"	44-54
Large	22-24.5"	54-60
X-Large	24.5-27.5"	60-67

18022 Cowan Street, Suite 200F, Irvine, CA 92614  
800.652.1135 Fax 800.652.1126 www.faststartortho.com

EXO is a trademark of VisionQuest Industries, Inc.  
FastStart is a registered trademark and service mark of VisionQuest Industries, Inc. © 2004

**FastStart**<sup>SM</sup>  
Orthopedic Products