



Post Operative Brace



Designed for immediate post operative use through complete range-of-motion.

- Maximum comfort and compliance
- Utilizes composite materials, which drastically reduces weight
- Precise polycentric hinge is easy to adjust to the desired range of motion
- Extension stops at 0, 5, 15, 30, 60, 90 degrees
- Flexion stops at 0, 15, 30, 60, 90 degrees
- Convenient snap-on hinge covers for easy access
- Slotted contoured arms hold strapping in place and conform to leg shape


FastStartSM
Orthopedic Products

800.652.1135



Post Operative Brace

Product Profile



ComfortLite Long

ComfortLite Mid/Long

ComfortLite Mid



Ankle attachment

These braces are universal size with a choice of 3 brace lengths: Mid = 18" or 44 cms. Mid/Long = 1.5" or 52.5 cms. Long = 24.5" or 60 cms.			
	Mid	Mid/Long	Long
Complete	07-30101-00-02	07-30101-00-03	07-30101-00-04
Accessories			
Ankle Attachment. Sizing (U.S. shoe size): Sm (M) 5-7 (F) 6 or less. Med (M) 7.5-10 (F) 6.5-10. Lg (M) 11-13 (F) 10+			
	Small	Medium	Large
Right	07-30104-00-22	07-30104-00-23	07-30104-00-24
Left	07-30104-00-12	07-30104-00-13	07-30104-00-14

18022 Cowan Street, Suite 200F, Irvine, CA 92614
 800.652.1135 Fax 800.652.1126 www.faststartortho.com

EXO is a trademark of VisionQuest Industries, Inc.
 FastStart is a registered trademark and service mark of VisionQuest Industries, Inc. © 2004

FastStartSM
 Orthopedic Products