

STx - Saunders Lumbar Traction Device™

Effective for either clinical or home use



Although the STx was designed for patient home use, its unique features are causing healthcare providers to use it as a clinical device!

What makes it so effective? Force is applied via a hand held pump, actively separating the lower half of the treatment surface. The active, friction-free force is highly efficient. Less pressure is transmitted through the pelvic belts, concentrating the traction force just where it's needed.

Features

- Friction-free track and actively moving treatment surface allow smooth and highly efficient application of lumbar traction
- Patented pump incorporates an easy "click into place" mechanism for pumping, sustaining and releasing the traction. It's easy to use and eliminates air leaks
- Gauge is calibrated in pounds and kilograms to provide feedback and promote consistency in treatment
- Weighs less than 30 pounds and can be set up in 30 seconds
- Storage compartment holds belts, pump, instructional video and manual
- Applies up to 200 pounds of traction
- Patient is held secure by comfortable thoracic and pelvic belts; most people can position themselves in the device without any help
- Safety feature: our unique pressure valve limits the amount of force to 200 pounds, and the traction force can be released easily and immediately by the patient if necessary

A physician's prescription is required for the Lumbar STx.

18022 Cowan Street, Suite 200F, Irvine, CA 92614
Fax 800.652.1126 www.faststartortho.com

800.652.1135

FastStart™
Orthopedic Products

FastStart is a registered trademark and service mark of VisionQuest Industries, Inc. © 2004
STx-Saunders Lumbar Traction Device is a trademark of The Saunders Group Physical Therapy Equipment

FS14032200