

## VQ™ Upper Limb Brace

Easy-to-use upper limb brace delivers superior patient comfort and protection.

- Range-of-motion and fixed position
- Available in both neutral and supination/pronation models to meet specific patient needs
- Optional over-the-shoulder humeral attachment for increased stability and protection



### PRODUCT PART NUMBERS

#### NEUTRAL

Part #	Size	Forearm Circumference	Bicep Circumference
01-30300-00-02	S	7.0-11.0 IN. (18-28 CM)	10.0-13.0 IN. (25-33 CM)
01-30300-00-03	M	10.0-13.0 IN. (25-33 CM)	12.0-15.0 IN. (30-38 CM)
01-30300-00-04	L	12.0-14.0 IN. (30-35 CM)	14.0-17.0 IN. (35-43 CM)
01-30300-00-05	XL	13.0-15.0 IN. (33-38 CM)	16.0- 19.0 IN. (41-48 CM)

#### SUPINATION / PRONATION

Part #	Size	Forearm Circumference	Bicep Circumference
07-30312-02	S	7.0-11.0 IN. (18-28 CM)	10.0-13.0 IN. (25-33 CM)
07-30312-03	M	10.0-13.0 IN. (25-33 CM)	12.0-15.0 IN. (30-38 CM)
07-30312-04	L	12.0-14.0 IN. (30-35 CM)	14.0-17.0 IN. (35-43 CM)
07-30312-05	XL	13.0-15.0 IN. (33-38 CM)	16.0- 19.0 IN. (41-48 CM)

#### HUMERAL ATTACHMENT (OVER THE SHOULDER)

Part #	Size
07-30301-02	S
07-30301-03	M
07-30301-04	L

TO ORDER CALL 800.652.1135



**INDICATIONS**

- Humeral and forearm fractures
- Post-operative shoulder management
- Biceps repair
- Ulnar nerve transposition

**CAUTION:** Applicable federal, state, and/or local regulations may restrict this device to sale by or on the order of a practitioner licensed by law of the State in which he/she practices to use or order the use of this device.



18011 Mitchell South, Irvine, CA 92614  
PHONE 800.266.6969 FAX 800.652.1126  
www.vqorthocare.com  
customersupport@vqorthocare.com